

ABL

Advanced
Biological
Laboratories



Improving Disease Management

DEEPCHEK®-HIV ASSAYS & SOFTWARE

End-to end Solutions for HIV Genotyping for Drug Resistance Assessment (SANGER & NGS)

PRESENTATION



Confidential & Proprietary ABL
Do Not Distribute – Strictly Private and Confidential

Presentation of ABL






Advanced Biological Laboratories (ABL), S.A., is a diagnostic and medical software company founded in 2000 as a spin-off from LIH (Luxembourg)


ABL's mission is to develop, manufacture and commercialize proprietary end-to-end diagnostic solutions for infectious diseases and oncology, in recurrent chronic diseases

The group is a 100% privately owned company, with offices in Luxembourg, France, Spain, South Africa and United States





Presentation of ABL Genotyping solutions








	SANGER	NGS
Turnaround time	~17h	~3-4 days
Hands-on-time	~40min	~1-3 h















HIV

HBV

HCV

TB




HPV

CMV

HSV

16s RNA




HIV Genotyping solutions



DEEPCHEK®-HIV ASSAYS & SOFTWARE

Workflow Overview



DEEPCHEK®- HIV ASSAYS

Highlights



- LoD: 1000 copies/mL (down to ~50-100 copied/mL)
- High analytical reproducibility and repeatability were evidenced by Percent Agreement being 100%
- Highly specific –pure & CRF strains
- From sample to drug resistance in ~17 hours for SANGER & 3-4 daya for NGS for 24 - 96 samples

DeepChek® Assay	Downstream sequencing instrument used with DeepChek® Assay	Device 2 used for agreement concordance	Nb. of samples tested	Concordance (%)
PR/RT + INT	Illumina MiSeq	Abbott® Dx - ViroSeq® HIV-1 Genotyping PR/RT + INT (Sanger)	23	100%
PR/RT + INT	Illumina MiSeq	Vela Dx - Sentosa® HIV-1 Genotyping PR/RT/INT (NGS)	18	100%



100% agreement with annual ANRS
quality control programs



100% agreement with annual QCMD
panel (RT+PR & INT) studies




DEEPCHEK®-HIV ASSAYS & SOFTWARE

ABL Plan for IVDR

HIV

Timelines to meet the
EU MDR and IVDR
compliance requirements




- Production of the IVD kits is already ISO 9001 certified
- Full design, development, production and support of IVD medical devices process is compliant with the ISO 13485:2016
- Since 2019, the Quality Management System (QMS) is developed to be meet the essential requirements of the new EU regulation for IVD medical devices (2017/746) for Class C and D products (highest classification)
- The Compliance team is coached and supported externally by a german group of experts in medical devices (reagents, instruments and software) to be ready for the certification
- Contacts with notified bodies to review and certify each IVD product according to the new EU regulation (ISO 13485 and technical file)
- ABL has a regulatory strategy for each product to be compliant with the new EU regulation, including the transition period



SANGER



SUBTYPING		COVERED POSITIONS	
PROT	B (94.95%)	RT	B (95.83%)
INT	B (95.53%)		
PROT 1-99		RT 1-327	
INT 32-289			
ANTIRETROVIRAL DRUG		INTERPRETATION	DEEPCHEK® HIV (Stanford 8.8)
Generic name	Assessment	Resistance mutations >20.00%	Resistance mutations between >3% and <20.00%
PI	Atazanavir/r	R	I54L (78.56%), G73S (76.1%), L90M (90.98%)
	Darunavir/r	I	I54L (78.56%)
	Fosamprenavir/r	R	I54L (78.56%), G73S (76.1%), L90M (90.98%)
	Indinavir/r	R	I54L (78.56%), G73S (76.1%), L90M (90.98%)
	Lopinavir/r	I	I54L (78.56%), G73S (76.1%), L90M (90.98%)
	Nelfinavir	R	I54L (78.56%), G73S (76.1%), L90M (90.98%)
	Saquinavir/r	R	I54L (78.56%), G73S (76.1%), L90M (90.98%)
NRTI	Tipranavir/r	S	I54L (78.56%)
	Abacavir	I	D67N (96.96%), M184V (99.35%)
	Didanosine	I	D67N (96.96%), M184V (99.35%)
	Emtricitabine	R	M184V (99.35%)
	Lamivudine	R	M184V (99.35%)
	Stavudine	S	D67N (96.96%), M184V (99.35%)
NNRTI	Tenofovir	S	D67N (96.96%), M184V (99.35%)
	Zidovudine	S	D67N (96.96%), M184V (99.35%)
	Doravirine	S	
	Elvitegravir	S	
	Etravirine	S	
INI	Nevirapine	S	
	Rilpivirine	S	
	Bictegravir	I	G140S (99.66%), Q148H (99.33%)
	Dolutegravir	I	G140S (99.66%), Q148H (99.33%)
S	Elvitegravir	R	G140S (99.66%), Q148H (99.33%)
	Raltegravir	R	G140S (99.66%), Q148H (99.33%)
Susceptible (S)		Potential low-level resistance (PLLR)	Low-level resistance (LLR)
		Intermediate resistant (IR)	High-level resistance (HLR)



VIROSCORE® & DEEPCHEK®- HIV SOFTWARE Highlights



Technology	Security	Main features	Additional features	Services
<ul style="list-style-type: none"> •Web-based (browser only) •Software & database •Local or Cloud (+HDS) hosting •Unlimited user accounts •Unlimited analyses 	<ul style="list-style-type: none"> •Data access restriction (pools, read-only mode...) •Logging of user accesses •Encrypted database •Reports validation 	<ul style="list-style-type: none"> •SANGER data management (AB1, FASTA) with embedded chromatogram editor •NGS data management (FASTQ) with dedicated pipeline •Genotyping (per sample & cumulative) •Subtyping •Virtual-phenotyping •Drug resistance through up-to-date guidelines •GSS determination and regiment ranking •Tropism •Reporting & labelling 	<ul style="list-style-type: none"> •Report customization •Contamination check •Quality control •Export (reports, FASTA, XML...) •Batch mode analysis •Data mining 	<ul style="list-style-type: none"> •Constant updates •Annual upgrades (versions) •Historical data import •LIMS integration •HIS integration •Support •Trainings

A WEB-BASED SYSTEM

The screenshot shows the ViroScore Premium Demo login interface. At the top, the ViroScore logo is displayed with the version number 3.28. Below the logo, it states 'VS Sanger Genotyping software and database solution' and 'VIROSCORE PREMIUM DEMO'. On the left, there is a support contact section with links to 'submit a support ticket', 'email customer support', and a support line number. On the right, there is a login form with fields for 'Login' (containing 'admin') and 'Password' (containing '*****'). Below these fields are radio buttons for 'User login' and 'Administrator login', with 'Administrator login' selected. There is a 'Sign-in' button and a link for 'Lost your password?'. At the bottom, a footer contains copyright information for 2018 and links to 'Terms of Use' and 'Privacy and Security Policy'.

VIROSCORE® & DEEPCHEK®- HIV SOFTWARE Highlights



Technology	Security	Main features	Additional features	Services
<ul style="list-style-type: none"> •Web-based (browser only) •Software & database ▶•Local or Cloud (+HDS) hosting •Unlimited user accounts •Unlimited analyses 	<ul style="list-style-type: none"> •Data access restriction (pools, read-only mode...) •Logging of user accesses •Encrypted database •Reports validation 	<ul style="list-style-type: none"> •SANGER data management (AB1, FASTA) with embedded chromatogram editor •NGS data management (FASTQ) with dedicated pipeline •Genotyping (per sample & cumulative) •Subtyping •Virtual-phenotyping •Drug resistance through up-to-date guidelines •GSS determination and regiment ranking •Tropism •Reporting & labelling 	<ul style="list-style-type: none"> •Report customization •Contamination check •Quality control •Export (reports, FASTA, XML...) •Batch mode analysis •Data mining 	<ul style="list-style-type: none"> •Constant updates •Annual upgrades (versions) •Historical data import •LIMS integration •HIS integration •Support •Trainings


SYSTEM AVAILABILITY




Cloud



Local Virtual Private Server



Secured Healthcare environment (HDS)



Local Hardware

VIROSCORE® & DEEPCHEK®- HIV SOFTWARE Highlights



Technology	Security	Main features	Additional features	Services
<ul style="list-style-type: none"> •Web-based (browser only) •Software & database •Local or Cloud (+HDS) hosting •Unlimited user accounts •Unlimited analyses 	<ul style="list-style-type: none"> •Data access restriction (pools, read-only mode...) •Logging of user accesses •Encrypted database •Reports validation 	<ul style="list-style-type: none"> •SANGER data management (AB1, FASTA) with embedded chromatogram editor •NGS data management (FASTQ) with dedicated pipeline •Genotyping (per sample & cumulative) •Subtyping •Virtual-phenotyping •Drug resistance through up-to-date guidelines •GSS determination and regment ranking •Tropism •Reporting & labelling 	<ul style="list-style-type: none"> •Report customization •Contamination check •Quality control •Export (reports, FASTA, XML...) •Batch mode analysis •Data mining 	<ul style="list-style-type: none"> •Constant updates •Annual upgrades (versions) •Historical data import •LIMS integration •HIS integration •Support •Trainings

SECURED DATA ACCESS

MULTI-SITE
COLLABORATIVE ACCESS
POSSIBLE


User Details			
First name *	<input type="text" value="Dimitri"/>	Last name *	<input type="text" value="Gonzalez"/>
Login *	<input type="text" value="dimitri"/>	Display name *	<input type="text" value="Dimitri (Demo)"/>
Email *	<input type="text" value="dimitri.gonzalez@ablsa.com"/>	Language	<input type="text" value="English"/>
Password *	<input type="password"/>	Other id	<input type="text"/>
Deactivated	<input type="checkbox"/>	<input type="button" value="Reset"/>	
User Privileges			
Role	<input type="text" value="Default"/>	Right *	<input type="text" value="show / add / edit / delete"/>
Data Access	<input type="text" value="All the patients"/>		
Administrator	<input checked="" type="checkbox"/>	Anonymizing Patient Data	<input type="checkbox"/>

VIROSCORE® & DEEPCHEK®- HIV SOFTWARE Highlights



Technology	Security	Main features	Additional features	Services
<ul style="list-style-type: none"> •Web-based (browser only) •Software & database •Local or Cloud (+HDS) hosting •Unlimited user accounts •Unlimited analyses 	<ul style="list-style-type: none"> •Data access restriction (pools, read-only mode...) •Logging of user accesses •Encrypted database •Reports validation 	<ul style="list-style-type: none"> •SANGER data management (AB1, FASTA) with embedded chromatogram editor •NGS data management (FASTQ) with dedicated pipeline •Genotyping (per sample & cumulative) •Subtyping •Virtual-phenotyping •Drug resistance through up-to-date guidelines •GSS determination and regimen ranking •Tropism •Reporting & labelling 	<ul style="list-style-type: none"> •Report customization •Contamination check •Quality control •Export (reports, FASTA, XML...) •Batch mode analysis •Data mining 	<ul style="list-style-type: none"> •Constant updates •Annual upgrades (versions) •Historical data import •LIMS integration •HIS integration •Support •Trainings

EMBEDDED CHROMATOGRAM EDITOR



VIROSCORE® & DEEPCHEK®- HIV SOFTWARE Highlights

HIV

Technology

- Web-based (browser only)
- Software & database
- Local or Cloud (+HDS) hosting
- Unlimited user accounts
- Unlimited analyses

Security

- Data access restriction (pools, read-only mode...)
- Logging of user accesses
- Encrypted database
- Reports validation

Main features

- SANGER data management (AB1, FASTA) with embedded chromatogram editor
- NGS data management (FASTQ) with dedicated pipeline
- Genotyping (per sample & cumulative)
- Subtyping
- Virtual-phenotyping
- Drug resistance through up-to-date guidelines
- GSS determination and regiment ranking
- Tropism
- Reporting & labelling

Additional features

- Report customization
- Contamination check
- Quality control
- Export (reports, FASTA, XML...)
- Batch mode analysis
- Data mining

Services

- Constant updates
- Annual upgrades (versions)
- Historical data import
- LIMS integration
- HIS integration
- Support
- Trainings

SINGLE or CUMULATIVE GENOTYPING

2017-09-26 2017-10-02 2017-10-13 RAPPORT CUMULATIF & HISTORIQUE - Viroscore+PLUS

Imprimer (html) Imprimer (couleur) Rapport PDF

CUMULATIVE & HISTORY REPORT - VIROSCORE PLUS (Mutations based)


Information sur l'échantillon		Analyse effectuée RT-PROT		Type de rapport	
ABLITE ID	Nom du patient	Nom du clinicien	Technique de séquençage	Date de rapport	Version du rapport
	Votre numéro d'échantillon		Type d'échantillon	05/02/2018 12:11:21	11.3
	Votre numéro de patient your_patient_84		Date de l'échantillon		
			Date du séquençage		

Mutations Dectectées (Séquence de référence HXB2) Mutations d'intérêt, en gras, basées sur Stanford v8.4 (mutation score = 0)

Mutation	Score
V81, V81V, V211, V211V, E28K, K32E, K32T, V35K, T39D, E53D, V60I, K65R , D67G, S66G, T69W, K70R , L74I, V75I , L100I, V106M , Y115F, I132L, I167V, K173T, M184V , Q196M, T215Y , K219E , P225H, L228R, T253S, T253T, V261Q, V261V, K275Q, R277K, Q278H	0% < x ≤ 25%
V35T, E36E, T39E, M41L , S48T, K103N , D123S, K173A, Q174K, I178M, R211K, V245Q, E291D, V292I, I293V	25% < x ≤ 50%
E36A, D177E, T200A, Q207E, L214F	50% < x ≤ 75%
T12S, I15V, L19T, T26S, T26T, E34Q, E35D, E35E, M46I , I54V, I62V, L63A, L63P, L76V , V82I, V82A , L89L, L89V	75% < x ≤ 100%
V31, L17H, K20R, R41K, M36I, S37N, H69K, L89M, I93L	0% < x ≤ 25%

Interprétation génotypiques (ANRS-27)

Classe de médicament	Médicament	ANRS 27
INRT	Zidovudine	R
	Didanosine	R
	Stavudine	R
	Lamivudine	R
	Emtricitabine	R
	Abacavir	R



VIROSCORE® & DEEPCHEK®- HIV SOFTWARE Highlights

HIV

Technology	Security	Main features	Additional features	Services
<ul style="list-style-type: none"> •Web-based (browser only) •Software & database •Local or Cloud (+HDS) hosting •Unlimited user accounts •Unlimited analyses 	<ul style="list-style-type: none"> •Data access restriction (pools, read-only mode...) •Logging of user accesses •Encrypted database •Reports validation 	<ul style="list-style-type: none"> •SANGER data management (AB1, FASTA) with embedded chromatogram editor •NGS data management (FASTQ) with dedicated pipeline •Genotyping (per sample & cumulative) •Subtyping •Virtual-phenotyping •Drug resistance through up-to-date guidelines •GSS determination and regimen ranking •Tropism •Reporting & labelling 	<ul style="list-style-type: none"> •Report customization •Contamination check •Quality control •Export (reports, FASTA, XML...) •Batch mode analysis •Data mining 	<ul style="list-style-type: none"> •Constant updates •Annual upgrades (versions) •Historical data import •LIMS integration •HIS integration •Support •Trainings

SUBTYPING


COMET

COMET HIV-1

(Context-based Modeling for Expeditious Typing)



PHYLOGENETIC ANALYSIS




VIROSCORE® & DEEPCHEK®- HIV SOFTWARE Highlights



Technology	Security	Main features	Additional features	Services
<ul style="list-style-type: none"> •Web-based (browser only) •Software & database •Local or Cloud (+HDS) hosting •Unlimited user accounts •Unlimited analyses 	<ul style="list-style-type: none"> •Data access restriction (pools, read-only mode...) •Logging of user accesses •Encrypted database •Reports validation 	<ul style="list-style-type: none"> •SANGER data management (AB1, FASTA) with embedded chromatogram editor •NGS data management (FASTQ) with dedicated pipeline •Genotyping (per sample & cumulative) •Subtyping ▶•Virtual-phenotyping •Drug resistance through up-to-date guidelines •GSS determination and regimen ranking •Tropism •Reporting & labelling 	<ul style="list-style-type: none"> •Report customization •Contamination check •Quality control •Export (reports, FASTA, XML...) •Batch mode analysis •Data mining 	<ul style="list-style-type: none"> •Constant updates •Annual upgrades (versions) •Historical data import •LIMS integration •HIS integration •Support •Trainings

VIRTUAL-PHENOTYPING



VIROSCORE® & DEEPCHEK®- HIV SOFTWARE Highlights



Technology	Security	Main features	Additional features	Services
<ul style="list-style-type: none"> •Web-based (browser only) •Software & database •Local or Cloud (+HDS) hosting •Unlimited user accounts •Unlimited analyses 	<ul style="list-style-type: none"> •Data access restriction (pools, read-only mode...) •Logging of user accesses •Encrypted database •Reports validation 	<ul style="list-style-type: none"> •SANGER data management (AB1, FASTA) with embedded chromatogram editor •NGS data management (FASTQ) with dedicated pipeline •Genotyping (per sample & cumulative) •Subtyping •Virtual-phenotyping •Drug resistance through up-to-date guidelines •GSS determination and regimen ranking •Tropism •Reporting & labelling 	<ul style="list-style-type: none"> •Report customization •Contamination check •Quality control •Export (reports, FASTA, XML...) •Batch mode analysis •Data mining 	<ul style="list-style-type: none"> •Constant updates •Annual upgrades (versions) •Historical data import •LIMS integration •HIS integration •Support •Trainings

DRUG RESISTANCE ASSESSMENT

- HIVDB/STANFORD
- ANRS
- REGA
- GENO2PHENO
- HIV-GRADE
- RIS
- RENAGENO

SANGER

Drug		ANRS 31 12-2020		STAN 9.0 03-2021
Zidovudine	R	D67N, L210W, T215Y, K219Q	R	D67N, L210W, T215Y, K219Q
Didanosine	Not available		R	D67N, T69D, L210W, T215Y, K219Q
Stavudine	Not available		R	D67N, T69D, L210W, T215Y, K219Q
Lamivudine	S		S	L210W, T215Y
Emtricitabine	S		S	L210W, T215Y
Abacavir	R	D67N, L210W, T215Y	I	D67N, L210W, T215Y, K219Q
Tenofovir	R	E44D, D67N, T69D, L210W, T215Y	I	D67N, L210W, T215Y, K219Q

NGS

	Algorithm	Sanger based sequencing	15.00%	10.00%	5.00%	3.00%
Abacavir	ANRS	I	I	I	I	I
	Grade	I	I	I	I	I
	Rega institute	S	S	S	S	S
	Stanford	LLR	LLR	LLR	LLR	LLR
Didanosine	ANRS	NA	NA	NA	NA	NA
	Grade	I	I	I	I	I
	Rega institute	S	S	S	S	S
	Stanford	LLR	LLR	LLR	LLR	LLR
Emtricitabine	ANRS	R	R	R	R	R
	Grade	R	R	R	R	R
	Rega institute	R	R	R	R	R
	Stanford	HLR	HLR	HLR	HLR	HLR

VIROSCORE® & DEEPCHEK®- HIV SOFTWARE Highlights



Technology	Security	Main features	Additional features	Services
<ul style="list-style-type: none"> •Web-based (browser only) •Software & database •Local or Cloud (+HDS) hosting •Unlimited user accounts •Unlimited analyses 	<ul style="list-style-type: none"> •Data access restriction (pools, read-only mode...) •Logging of user accesses •Encrypted database •Reports validation 	<ul style="list-style-type: none"> •SANGER data management (AB1, FASTA) with embedded chromatogram editor •NGS data management (FASTQ) with dedicated pipeline •Genotyping (per sample & cumulative) •Subtyping •Virtual-phenotyping •Drug resistance through up-to-date guidelines •GSS determination and regimen ranking •Tropism •Reporting & labelling 	<ul style="list-style-type: none"> •Report customization •Contamination check •Quality control •Export (reports, FASTA, XML...) •Batch mode analysis •Data mining 	<ul style="list-style-type: none"> •Constant updates •Annual upgrades (versions) •Historical data import •LIMS integration •HIS integration •Support •Trainings

GSS & REGIMEN RANKING

Recommendation	Combination Regimen	Algorithm	Sanger based sequencing	15.00% ↓	10.00%	5.00%	3.00%
AI	Emtricitabine Tenofovir Dolutegravir	ANRS	1 / 3	1 / 3 *	1 / 3 *	1 / 3 *	1 / 3 *
		Grade	0.5 / 3	0.5 / 3 *	0.5 / 3 *	0.5 / 3 *	0.5 / 3 *
		Rega institute	1 / 3	2 / 3	2 / 3	2 / 3	2 / 3
		Stanford ↓	1 / 3	2 / 3	2 / 3	2 / 3	2 / 3
AI	Emtricitabine Tenofovir Atazanavir/r	ANRS	2 / 3	2 / 3	2 / 3	2 / 3	2 / 3
		Grade	1.5 / 3	1.5 / 3	1.5 / 3	1.5 / 3	1.5 / 3
		Rega institute	2 / 3	2 / 3	2 / 3	2 / 3	2 / 3
		Stanford ↓	2 / 3	2 / 3	2 / 3	2 / 3	2 / 3
AI	Lamivudine Tenofovir Atazanavir/r	ANRS	2 / 3	2 / 3	2 / 3	2 / 3	2 / 3
		Grade	1.5 / 3	1.5 / 3	1.5 / 3	1.5 / 3	1.5 / 3
		Rega institute	2 / 3	2 / 3	2 / 3	2 / 3	2 / 3
		Stanford ↓	2 / 3	1 / 3	2 / 3	2 / 3	2 / 3

VIROSCORE® & DEEPCHEK®- HIV

SOFTWARE Highlights



Technology	Security	Main features	Additional features	Services
<ul style="list-style-type: none"> •Web-based (browser only) •Software & database •Local or Cloud (+HDS) hosting •Unlimited user accounts •Unlimited analyses 	<ul style="list-style-type: none"> •Data access restriction (pools, read-only mode...) •Logging of user accesses •Encrypted database •Reports validation 	<ul style="list-style-type: none"> •SANGER data management (AB1, FASTA) with embedded chromatogram editor •NGS data management (FASTQ) with dedicated pipeline •Genotyping (per sample & cumulative) •Subtyping •Virtual-phenotyping •Drug resistance through up-to-date guidelines •GSS determination and regiment ranking •Tropism •Reporting & labelling 	<ul style="list-style-type: none"> •Report customization •Contamination check •Quality control •Export (reports, FASTA, XML...) •Batch mode analysis •Data mining 	<ul style="list-style-type: none"> •Constant updates •Annual upgrades (versions) •Historical data import •LIMS integration •HIS integration •Support •Trainings

TROPISM DETERMINATION

DeepChek® - HIV Co-receptor tropism

X4

X4

99.909363818714%

R5

0.090636181285916%

Activity of CCR5 antagonist anticipated ?

YES

NO

☐
☒

R5

Virus uses CCR5 co-receptors to enter CD4+ cells

D/M

Virus uses CCR5 and/or CXCR4 co-receptors to enter CD4+ cells

X4

Virus uses CXCR4 co-receptors to enter CD4+ cells

VIROSCORE® & DEEPCHEK®- HIV SOFTWARE Highlights



Technology	Security	Main features	Additional features	Services
<ul style="list-style-type: none"> •Web-based (browser only) •Software & database •Local or Cloud (+HDS) hosting •Unlimited user accounts •Unlimited analyses 	<ul style="list-style-type: none"> •Data access restriction (pools, read-only mode...) •Logging of user accesses •Encrypted database •Reports validation 	<ul style="list-style-type: none"> •SANGER data management (AB1, FASTA) with embedded chromatogram editor •NGS data management (FASTQ) with dedicated pipeline •Genotyping (per sample & cumulative) •Subtyping •Virtual-phenotyping •Drug resistance through up-to-date guidelines •GSS determination and regimen ranking •Tropism •Reporting & labelling 	<ul style="list-style-type: none"> •Report customization •Contamination check •Quality control •Export (reports, FASTA, XML...) •Batch mode analysis •Data mining 	<ul style="list-style-type: none"> •Constant updates •Annual upgrades (versions) •Historical data import •LIMS integration •HIS integration •Support •Trainings

REPORT VALIDATION



Patient/Sample information

ABL/TE ID

Patient name

Patient birth date

Patient viral load Missing_data

Your sample ID Test demo 2

Your patient ID

Analysis RT-PROT-INT

Clinician name Dr toto

Genotyping method DeepChek® SingleRound

Type of sample Plasma

Date of sample 20/08/2020

Date of sequencing 21/08/2020

Type of report Initial

Data report edition 06/07/2021,10:29:19

Data report version 12.4

Report label

Validated - (TeSuperAdmin 06/07/2021 10:30:09)

Subtype:

Reverse transcriptase B (98 % similarity)

Protease B (79 % similarity)

Integrase B (100 % similarity)

All Mutations Detected (HXB2 reference Sequence) Resistance mutations in bold based on Stanford 9.0

Reverse transcriptase	K20R, V35I, E44D, D67N, T69D , V90I, V106I , V118I, I135T, A158S, S162C, T165I, Y181C, G190A, L210W, T215Y , D218E, K219Q , L228H, V245E, E248D, P272A, R277K, T286A, A288T, I293V, E297R
Protease	V3I, L10V, V11L , I13V, K14R, I15V, K20T, V32I, L33F , E35D, M36I, S37N, R41K, I54L, Q58E , D60E, Q61E, I62V, L63P, A71V, G73N, I84V, L89V, L90M
Integrase	V31L, del33, del34, del35, del36, del37, V72I, G123S, A124T, T125V, R127K, N232D, L234I, D253E, D256E

VIROSCORE® & DEEPCHEK®- HIV SOFTWARE Highlights



Technology	Security	Main features	Additional features	Services
<ul style="list-style-type: none"> •Web-based (browser only) •Software & database •Local or Cloud (+HDS) hosting •Unlimited user accounts •Unlimited analyses 	<ul style="list-style-type: none"> •Data access restriction (pools, read-only mode...) •Logging of user accesses •Encrypted database •Reports validation 	<ul style="list-style-type: none"> •SANGER data management (AB1, FASTA) with embedded chromatogram editor •NGS data management (FASTQ) with dedicated pipeline •Genotyping (per sample & cumulative) •Subtyping •Virtual-phenotyping •Drug resistance through up-to-date guidelines •GSS determination and regiment ranking •Tropism •Reporting & labelling 	<ul style="list-style-type: none"> •Report customization •Contamination check •Quality control •Export (reports, FASTA, XML...) •Batch mode analysis •Data mining 	<ul style="list-style-type: none"> •Constant updates •Annual upgrades (versions) •Historical data import •LIMS integration •HIS integration •Support •Trainings

CONTAMINATION CHECK



Screening for similar strains



Region	Similar strains			
Reverse transcriptase	No significant homology detected compared to the strains stored in the database.			
Protease	No significant homology detected compared to the strains stored in the database.			
	Sample ID	Sample date	Date performed	Identities
	1	10/02/2021	04/03/2021 13:28:31	96.40%
	2	26/09/2019	24/10/2019 10:56:31	96.15%
	3	19/05/2021	07/06/2021 15:41:50	95.87%
	4	26/03/2020	19/10/2020 13:24:58	95.76%
	5	22/07/2020	16/10/2020 13:44:56	95.88%
Integrase	6	03/01/2020	27/01/2020 11:18:17	96.10%
	7	08/02/2021	15/03/2021 11:17:11	96.24%
	8	06/02/2020	03/03/2020 16:33:41	95.62%
	9	27/05/2021	18/06/2021 11:56:32	95.76%
	10	18/05/2021	07/06/2021 11:45:27	95.61%
		14/04/2021	23/06/2021 10:56:38	95.23%

Only the 10 first matches scored with an identity > 95% are displayed.

VIROSCORE® & DEEPCHEK®- HIV SOFTWARE Highlights

HIV

Technology

- Web-based (browser only)
- Software & database
- Local or Cloud (+HDS) hosting
- Unlimited user accounts
- Unlimited analyses

Security

- Data access restriction (pools, read-only mode...)
- Logging of user accesses
- Encrypted database
- Reports validation

Main features

- SANGER data management (AB1, FASTA) with embedded chromatogram editor
- NGS data management (FASTQ) with dedicated pipeline
- Genotyping (per sample & cumulative)
- Subtyping
- Virtual-phenotyping
- Drug resistance through up-to-date guidelines
- GSS determination and regimen ranking
- Tropism
- Reporting & labelling

Additional features

- Report customization
- Contamination check
- Quality control
- Export (reports, FASTA, XML...)
- Batch mode analysis
- Data mining

Services

- Constant updates
- Annual upgrades (versions)
- Historical data import
- LIMS integration
- HIS integration
- Support
- Trainings


QUALITY CONTROL

SANGER

Sequence	Quality				
REVERSE TRANSCRIPTASE	Good				
PROTEASE	Good				
INTEGRASE	Good				
GP120/V3Loop					
GP41					

	RT	PROT	INT	GP120/V3Loop	GP41
Length (nt)	936	297	778		
	(Good)	(Good)	(Good)		
Stop codons	0	0	0		
Insertions	0	0	0		
Deletions	0	0	5		
Frame shift	0	0	1		
	(Good)	(Good)	(Good)		
Ambiguity (2 bases/nt)	0	0	0		
Ambiguity (3 bases/nt)	0	0	0		
	(Good)	(Good)	(Good)		
Ambiguity (4 bases/nt)	0	0	0		
Bad characters	0	0	0		
Minimal position	1	1	30		
Maximal position	312	99	289		
Ambiguity 2 / resistance codon	0	0	0		
Ambiguity 3 / resistance codon	0	0	0		
Ambiguity 4 / resistance codon	0	0	0		
Sequence ambiguity (%)	0.00	0.00	0.00		

NGS



VIROSCORE® & DEEPCHEK®- HIV SOFTWARE Highlights



Technology	Security	Main features	Additional features	Services
<ul style="list-style-type: none"> •Web-based (browser only) •Software & database •Local or Cloud (+HDS) hosting •Unlimited user accounts •Unlimited analyses 	<ul style="list-style-type: none"> •Data access restriction (pools, read-only mode...) •Logging of user accesses •Encrypted database •Reports validation 	<ul style="list-style-type: none"> •SANGER data management (AB1, FASTA) with embedded chromatogram editor •NGS data management (FASTQ) with dedicated pipeline •Genotyping (per sample & cumulative) •Subtyping •Virtual-phenotyping •Drug resistance through up-to-date guidelines •GSS determination and regiment ranking •Tropism •Reporting & labelling 	<ul style="list-style-type: none"> •Report customization •Contamination check •Quality control •Export (reports, FASTA, XML...) •Batch mode analysis •Data mining 	<ul style="list-style-type: none"> •Constant updates •Annual upgrades (versions) •Historical data import •LIMS integration •HIS integration •Support •Trainings

DATA EXPORT



- ☐ Name
- ☒ FASTA
- ☒ PDF
- ☒ VS_Database_structure.doc
- ☒ VS_HIV.csv
- ☒ VS_HIV_GSS.csv
- ☒ VS_HIV_has_projects.csv
- ☒ VS_HIVclinical.csv
- ☒ VS_HIVinterpretations.csv
- ☒ VS_HIVmut_comments.csv
- ☒ VS_HIVquality.csv
- ☒ VS_HIVreports.csv
- ☒ VS_HIVsubtype.csv
- ☒ VS_HIVtreated.csv
- ☒ VS_projects.csv


VIROSCORE® & DEEPCHEK®- HIV SOFTWARE Highlights



Technology	Security	Main features	Additional features	Services
<ul style="list-style-type: none"> •Web-based (browser only) •Software & database •Local or Cloud (+HDS) hosting •Unlimited user accounts •Unlimited analyses 	<ul style="list-style-type: none"> •Data access restriction (pools, read-only mode...) •Logging of user accesses •Encrypted database •Reports validation 	<ul style="list-style-type: none"> •SANGER data management (AB1, FASTA) with embedded chromatogram editor •NGS data management (FASTQ) with dedicated pipeline •Genotyping (per sample & cumulative) •Subtyping •Virtual-phenotyping •Drug resistance through up-to-date guidelines •GSS determination and regiment ranking •Tropism •Reporting & labelling 	<ul style="list-style-type: none"> •Report customization •Contamination check •Quality control •Export (reports, FASTA, XML...) •Batch mode analysis •Data mining 	<ul style="list-style-type: none"> •Constant updates •Annual upgrades (versions) •Historical data import •LIMS integration •HIS integration •Support •Trainings

BATCH ANALYSIS
LARGESCALE

SANGER



NGS

Configuration de l'analyse

Données d'entrée	RT	PROT	INT	GP41	GP120	Votre ID d'échantillon *	Date de l'échantillon *	Charge virale (copies/mL)
Raw Seq Sample2INT	X					Sample2INT	dd/mm/yyyy	
Raw Seq Sample2	X					Sample2	dd/mm/yyyy	
Raw Seq Sample1	X					Sample1	dd/mm/yyyy	


Select an item from your dataset & Drag here to analyse

VIROSCORE® & DEEPCHEK®- HIV SOFTWARE Highlights



Technology	Security	Main features	Additional features	Services
<ul style="list-style-type: none"> •Web-based (browser only) •Software & database •Local or Cloud (+HDS) hosting •Unlimited user accounts •Unlimited analyses 	<ul style="list-style-type: none"> •Data access restriction (pools, read-only mode...) •Logging of user accesses •Encrypted database •Reports validation 	<ul style="list-style-type: none"> •SANGER data management (AB1, FASTA) with embedded chromatogram editor •NGS data management (FASTQ) with dedicated pipeline •Genotyping (per sample & cumulative) •Subtyping •Virtual-phenotyping •Drug resistance through up-to-date guidelines •GSS determination and regiment ranking •Tropism •Reporting & labelling 	<ul style="list-style-type: none"> •Report customization •Contamination check •Quality control •Export (reports, FASTA, XML...) •Batch mode analysis •Data mining 	<ul style="list-style-type: none"> •Constant updates •Annual upgrades (versions) •Historical data import •LIMS integration •HIS integration •Support •Trainings


DATA MINING



VIROSCORE® & DEEPCHEK®- HIV SOFTWARE Highlights



Technology	Security	Main features	Additional features	Services
<ul style="list-style-type: none"> •Web-based (browser only) •Software & database •Local or Cloud (+HDS) hosting •Unlimited user accounts •Unlimited analyses 	<ul style="list-style-type: none"> •Data access restriction (pools, read-only mode...) •Logging of user accesses •Encrypted database •Reports validation 	<ul style="list-style-type: none"> •SANGER data management (AB1, FASTA) with embedded chromatogram editor •NGS data management (FASTQ) with dedicated pipeline •Genotyping (per sample & cumulative) •Subtyping •Virtual-phenotyping •Drug resistance through up-to-date guidelines •GSS determination and regimen ranking •Tropism •Reporting & labelling 	<ul style="list-style-type: none"> •Report customization •Contamination check •Quality control •Export (reports, FASTA, XML...) •Batch mode analysis •Data mining 	<ul style="list-style-type: none"> •Constant updates •Annual upgrades (versions) •Historical data import •LIMS integration •HIS integration •Support •Trainings





- Advanced Biological Laboratories S.A.
- 52- 54 Av du X Septembre, Luxembourg
- contact@ablsa.com
- <https://www.ablsa.com>
- TEL : (+352) 26 38 96 761
- FAX : (+352) 26 38 96 76 31

<http://www.ablsa.com>