


# UltraGene Assay SARS-CoV-2 Universalis (V2.X)

qPCR Assay to detect specifically the SARS-CoV-2 RNA including **current circulating lineages** without drop-out or lack of sensitivity



The UltraGene Assay SARS-CoV-2 Universalis (V2.X) is a real-time RT-PCR test intended for the qualitative detection of RNA from the SARS-CoV-2 in upper respiratory specimens from individuals suspected of COVID-19 by their healthcare provider

*Positive results are indicative of the presence of SARS-CoV-2 RNA; clinical correlation with patient history and other diagnostic information is necessary to determine patient infection status.*

## Workflow

**Sample Collection**

**Manual**   **Automated**

~25min   ~30min

**Extraction**

~20min   ~1h10

**PCR Set Up**   **Run**

**rt-qPCR**

**Data Analysis**

# Overview

Features	UltraGene Assay SARS-CoV-2 Universalis V2.X
Formats	<ul style="list-style-type: none"> <li>50 or 1000 tests per kit</li> <li>Multiplex format in 1 well per sample</li> </ul>
Content	<ul style="list-style-type: none"> <li>Master Mix 2X, RT Mix, Enzyme Mix and Magnesium Sulfate</li> <li>Specific Primers &amp; Probes for each targeted gene: N, ORF1ab (RdRp)</li> <li>Positive control</li> <li>Negative control</li> <li>Internal control</li> <li>PCR grade nuclease-free water</li> </ul>
Specimen	<ul style="list-style-type: none"> <li>Upper respiratory specimens (nasal, mid-turbinate, nasopharyngeal, oropharyngeal swab specimens and nasopharyngeal wash/aspirate or nasal aspirate and bronchoalveolar lavage (BAL) fluid specimens) from individuals suspected of COVID-19 by their healthcare provider</li> </ul>
Extraction	<ul style="list-style-type: none"> <li>Magnetic beads</li> <li>Validated for use with:                             <ul style="list-style-type: none"> <li>Seegene STARlet instrument (Cat. #67930-03) with Hamilton Run Control software (version 4.5.0.7977)</li> <li>Seegene STARMag 96 X 4 Universal Cartridge Kit (Cat. #744300.4.UC384)</li> <li>Any laboratory validated instrument for RNA extraction and purification using magnetic-bead technology shall work</li> </ul> </li> </ul>
Compatibility with qPCR platforms	<ul style="list-style-type: none"> <li>Any qPCR instrument compatible with the FAM, HEX, ROX and Cy5 channels.</li> <li>Validated for use with:                             <ul style="list-style-type: none"> <li>C1000 Dx Thermal Cycler (Bio-Rad, Catalog # 1841000-IVD with CFX96-Dx ORM software version 3.1) (C1000 CFX96),</li> <li>QuantStudio 5 Real-Time PCR Instrument (96-Well 0.1mL Block) (Applied Biosystems, Catalog #A28133, Design &amp; Analysis Software 1.5.2 / Firmware Version 1.4.0)</li> </ul> </li> </ul>
Hands-on-time	<ul style="list-style-type: none"> <li>Less than 20 minutes for up to 48 samples</li> <li>Less than 40 minutes for up to 96 samples</li> </ul>
Limit of detection (LOD)	<ul style="list-style-type: none"> <li>1.15 TCID<sub>50</sub>/mL for SARS-CoV-2 with the C1000 CFX96 Dx Thermal Cycler and the QuantStudio 5 Real-Time PCR Instruments.</li> </ul>
Clinical performance evaluation	<ul style="list-style-type: none"> <li>The positive (PPA) and negative (NPA) percent agreements between the UltraGene Assay SARS-CoV-2</li> <li>Universalis V2.X and CE-IVD EUA test are for the C1000 CFX96:                             <ul style="list-style-type: none"> <li>PPA = 28/28 = 100% (95% C.I. = 84.98% - 100%)</li> <li>NPA = 30/30 = 100% (95% C.I. = 85.86% - 100%)</li> </ul> </li> </ul>
Clinical Accuracy	<ul style="list-style-type: none"> <li>96%</li> </ul>



UltraGene Assay SARS-CoV-2 Universalis (V1.x) (50 tests)	169B50
UltraGene Assay SARS-CoV-2 Universalis (V1.x) (1000 tests)	169B1000



- Advanced Biological Laboratories S.A.
- 52-54 avenue X Septembre
- L-2550 Luxembourg, Luxembourg
- [contact@ablsa.com](mailto:contact@ablsa.com)
- <https://www.ablsa.com>
- TEL : (+352) 26 38 96 761
- FAX : (+352) 26 38 96 76 31