

of HIV-1 and Influenza A/B viruses Using Nanopore Technology



Rezak Drali¹, Amira Doudou¹, Chalom Sayada², Sofiane Mohamed¹
¹ABL Diagnostics, FRANCE; ²ABL Diagnostics SA, LUXEMBOURG



www.linkedin.com/company/abl-sa
Contact information: s.mohamed@ablsa.com

Introduction

The Human Immunodeficiency Virus-1 (HIV-1) and the Influenza Virus are RNA viruses responsible for serious illness, causing hundreds of thousands of deaths worldwide each year. The development of high-performance surveillance tools, such as whole genome sequencing, has made it possible to make significant advances in the knowledge of these pathogens with a view to controlling and/or eradicating them in the medium and long term. In this work, we show that simultaneous sequencing of the whole genomes of these two viruses using Nanopore technology offered an excellent alternative to expensive technologies.


Methods

Twelve samples from PCR-positive patients (5 for HIV-1 and 7 for influenza) were prepared with the QIAamp Mini (Qiagen) kit. OneStep RT-PCR was performed with the ABL Diagnostics DeepChek® Whole Genome HIV-1 and DeepChek® Whole Genome Influenza A/B kits, respectively. Whole genome sequencing with the MinION Nanopore device was performed after multiplexing the 12 samples included in this study. To evaluate the accuracy of Nanopore sequencing, a set of sequences was also subjected to Illumina iSeq100 sequencing. DeepChek® (ABL) software was used to assemble genomes, detect mutations, and build phylogenetic trees.

Results

The sequences generated by the MinION and Illumina technologies were comparatively similar and consistent with the reference genomes of the HIV-1 and Influenza A/B viruses analyzed here. Coverage of the sequenced genomes was complete for all 12 targeted samples, namely, the five HIV-1 samples (subtypes: two B, two C and one O2AG) and the seven influenza samples (two A/H1N1, two A/H3N2 and three B/Victoria). The identified subtype and mutations (>20%) found in Illumina and Nanopore were similar.

Results of multiplex Whole Genome Sequencing
Minion Oxford Nanopore vs Illumina



Conclusions

Nanopore technology, which enables multiplex sequencing of complete genomes of dangerous viruses such as HIV-1 and Influenza A/B, offers a real alternative to complex and expensive sequencing methods, especially for developing countries. In addition to being specific and compact, Nanopore technology is easy to use. This will inevitably have an impact on the monitoring of antiviral resistance as well as on epidemiological and clinical studies soon.

