

Optimizing Clinical NGS: Comparative Evaluation of Illumina iSeq100 and Complete Genomics G10-FR* Using DeepChek® Workflows


Hui Ren¹, David Du, Charles Bao¹, Sam Crawford¹, Rezak Drali², Jessica Arliaud², Margaux Mercadal², Aitor Modol³, Giuseppina Zuco², Chalom Sayada⁴, Ronan Boulmé², Dimitri Gonzalez², Laurent Deblir⁴, Matthieu Barralon⁴, Sofiane Mohamed²
¹ Complete Genomics, ²ABL Diagnostics, France ; ³ABL Spain ; ⁴ABL SA, Luxembourg

**The G10-FR is marketed as E25 Flash in some territories*

ID046

Background

Next-generation sequencing (NGS) is transforming infectious disease diagnostics by enabling comprehensive genotyping, resistance detection, and microbial identification^{1,2}. However, accessibility is often limited by cost, workflow complexity, and platform availability. While Illumina platforms are widely adopted³, newer systems like the Complete Genomics G10-FR* offer rapid sequencing with simplified workflows. In this study, we compared the performance of the Illumina iSeq100 and the Complete Genomics G10-FR (marketed as DNBSEQ-E25 Flash in some territories) using validated DeepChek® (ABL) kits and MicrobioChek® software, with a focus on pooling multiple applications in a single run to reduce per-sample costs.



Complete Genomics DNA Nanoball Sequencing:
☒ No clonal error
☒ No index hopping
☒ Low duplication rates

DeepChek® (ABL, Luxembourg) is a CE marked: The pipelines consist of 10 major steps which are :
1) Data cleaning, 2) Subtyping, 3) Tropism analysis, 4) Alignment, 5) Post-alignment cleaning, 6) Consensus creation, 7) Variant calling, 8) Expert system filtering, 9) Drug resistance calculation, 10) Reporting

- ☒ DeepChek® HIV version 2.0 / expert system 2.3 /algorithms 13.1
- ☒ ANRS 33 10-2022 / Grade 2021 9-2021 / HIVdb 9.4 12-2022
- ☒ Classification of mutations of interest ANRS 33
- ☒ Microbiochek identification

Figure 1. End-to end Solution for HIV-1 Genotyping and Drug Resistance for Routine Diagnostic Sequencing.

Results

Both platforms generated high-quality sequencing data with high Q-scores and full genome coverage across applications. The Complete Genomics G10-FR* produced a significantly higher number of reads per run (~25M) compared to the Illumina iSeq100 (~4M), enabling greater depth and broader detection capacity, particularly in pooled workflows.




Figure 2. Comparison of the number of HIV sequences generated by the sequencing platforms Complete Genomics G10-FR* versus iSeq100 (Illumina). 121A/122A : target specific kit. 170B: whole genome kit.

Phylogenetic analyses demonstrated strong concordance between sequences obtained on both platforms.




Figure 3. Phylogenetic trees constructed using consensus sequences of protease (a), reverse transcriptase (b), and integrase (c) targets obtained from sequencing on the NGS platforms Complete Genomics G10-FR* versus Illumina iSeq100.

HIV genotyping showed >99% agreement, with only two low-frequency mutations (<5%) showing discordance between platforms: M46V (3.32%) detected only by G10-FR* in the protease region, and D67N (3.25%) detected only by iSeq100 in the reverse transcriptase region. Importantly, these minor discrepancies (<5%) did not affect drug resistance interpretation thresholds. The complete workflow from DNA to data analysis takes approximately 24 hours on the Illumina iSeq100, including library preparation and sequencing. In comparison, Complete Genomics G10-FR* can deliver results in just 5-9 hours, offering a 50-70% faster turnaround.

Table 1. Bioinformatics analysis using MicrobioChek (ABL) to identify HIV-1 drug resistances obtained from Complete Genomics G10-FR* versus iSeq100 (Illumina) platforms.

N°	Sequencer	Subtype						Mutation						
		PR	%	RT	%	INT	%	PR	%	RT	%	INT	%	
QCMD2	MGI E25	95_02B	96.63	02_AG	96.15	B	96.66	NA	NA	NA	NA	NA	NA	
	iSeq 100	95_02B	96.63	02_AG	95.15	B	96.66	NA	NA	NA	NA	NA	NA	
QCMD1	G10-FR	48_01B	93.27	B	99.58	08_BC	97.09	43K = T 12.56 48M = I 95.95 54I = V 98.41 82V = A 95.76 90L = M 96.92 43K = T 12.07 46M = I 96.3 46M = V 3.32 54I = V 96.63 82V = A 91.35 90L = M 96.33	NA	NA	NA	NA		
	iSeq 100	48_01B	93.27	B	99.58	08_BC	97.09	NA	NA	NA	NA	NA	NA	
QCMD2	G10-FR	95_02B	96.63	02_AG	94.86	B	99.45	NA	NA	190G = E 4.35 190G = E 4.59 41M = L 99.35 67D = N 98.56 69T = D 99.31 98A = G 99.17 184M = I 98.6 188Y = L 98.05 190G = A 99.44 210I = W 99.05 215T = Y 98.29 41M = L 99.63 67D = N 99.01 69T = D 98.08 98A = G 99.18 184M = I 99.13 188Y = L 99.37 190G = A 99.48 210I = W 98.85 215T = Y 99.17	4.35	NA	NA	NA
	iSeq 100	95_02B	96.63	02_AG	94.86	B	99.45	NA	NA	NA	NA	NA	NA	
QCMD3	G10-FR	C	97.31	C	94.4	B	99.45	88N = S 10.65 98A = G 99.17 184M = I 98.6 188Y = L 98.05 190G = A 99.44 210I = W 99.05 215T = Y 98.29 41M = L 99.63 67D = N 99.01 69T = D 98.08 98A = G 99.18 184M = I 99.13 188Y = L 99.37 190G = A 99.48 210I = W 98.85 215T = Y 99.17	NA	NA	NA	NA		
	iSeq 100	C	97.31	C	94.4	B	99.45	88N = S 12.1 98A = G 99.18 184M = I 99.13 188Y = L 99.37 190G = A 99.48 210I = W 98.85 215T = Y 99.17	NA	NA	NA	NA		
QCMD4	G10-FR	146_BC	93.92	57_BC	94.93	B	99.58	46M = I 95.92 54I = V 98.41 82V = A 91.8 46M = I 96.01 54I = V 96.06 82V = A 90.55	185M = V 98.08 NA	143Y = H 4.98 NA	NA	4.68		
	iSeq 100	146_BC	93.92	57_BC	94.93	B	99.58	NA	NA	NA	NA	NA	NA	
QCMD5	G10-FR	B	99.66	B	96.14	121_0107	96.77	NA	NA	40E = F 96.43 41M = L 100 184M = V 96.88 210I = W 97.56 215T = Y 100 40E = F 100 41M = L 100 184M = V 100 210I = W 100 215T = Y 100	NA	NA	NA	NA
	iSeq 100	B	99.66	154_0755	98.21	121_0107	96.77	NA	NA	NA	NA	NA	NA	

Microbial identification via 16S/ITS was consistent across platforms at species/strain level.

Table 2. Bioinformatics analysis using MicrobioChek (ABL) to identify bacteria and fungi obtained from Complete Genomics G10-FR* platform.

N°ABL	PCR	Expected result	Result Complete Genomics G10-FR - MicrobioChek	Accession
60	16Sv3	Bordetella pertussis A639	Bordetella pertussis strain NCTC13251 chromosome 1	NZ_LR590467.1
62	16Sv3	Haemophilus influenzae type b (Eagan)	Haemophilus influenzae isolate KR21 chromosome 1, complete sequence	NZ_OV040584.1
64	16Sv3	Streptococcus pyogenes (2016)	Streptococcus pyogenes strain NCTC5318 chromosome 1	NZ_LR560683.1
67	16Sv3	Pseudomonas aeruginosa clinical isolate	Pseudomonas aeruginosa isolate MINF_3A-sc-2280432 chromosome 1, complete sequence	NZ_LR888867.1
69	16Sv3	Streptococcus salivarius (T127)	Streptococcus salivarius strain NCTC8618 chromosome 1, complete sequence	NZ_LR134274.1
5A	ITS	Staphylococcus aureus	Staphylococcus aureus strain 5720190863 isolate 5720190863 chromosome 1, complete sequence	NZ_OV344719.1
65	ITS - Primers 2	Saccharomyces cerevisiae recombinant	Saccharomyces cerevisiae S288C chromosome XII, complete sequence	NC_001444.5
66	ITS - Primers 2	Candida albicans 2006	Candida albicans SC5314 chromosome R, complete sequence	NC_032066.1
65	ITS - Primers 2	Saccharomyces cerevisiae S288C recombinant	Saccharomyces cerevisiae S288C chromosome XII, complete sequence	NC_001444.5
65	ITS - Primers 2	Candida albicans 2006	Candida albicans SC5314 chromosome R, complete sequence	NC_032066.1
63	16Sv3 (100ng)	Streptococcus pneumoniae (19F - 2022)	Streptococcus pneumoniae strain RMV7 isolate RMV7/dm1 chromosome 1, complete sequence	NZ_OV904788.1

Conclusion

DeepChek® NGS workflows are polyvalent, platform-independent, and compatible with a multiplexed, pooled sequencing approach, which significantly enhances cost-efficiency without sacrificing analytical performance. The Complete Genomics G10-FR* emerges as a robust alternative to traditional systems like Illumina iSeq100, making routine clinical NGS faster, more accessible and scalable for laboratories handling diverse diagnostic applications.

References

- Vellas, et al. Comparison of short-read and long-read next-generation sequencing technologies for determining HIV-1 drug resistance. J Med Virol. 2024
- Martinez, et al. Advancing Drug Resistance Detection: Comparative Analysis Using Short-Read and Long-Read NGS Technologies. LabMed 2025
- Mohamed, et al. From Capillary Electrophoresis to Deep Sequencing: An Improved HIV-1 Drug Resistance Assessment Solution Using In Vitro Diagnostic (IVD) Assays and Software. Viruses. 2023